

mRNA Data Analysis Pipeline

Quality Assessment

FastQC



Babraham Bioinformatics

Trimming

Skewer

Quality Assessment

FastQC



Babraham Bioinformatics

Mapping to a Reference

STAR



Integrative
Genomics
Viewer

Visualization

Counting reads per gene

HTSeq



OPEN SOURCE SOFTWARE FOR BIOINFORMATICS


Differential Gene Expression

DESeq2

Submit to Public Repository



What if you don't have a reference?



Quality Assessment

Quality Assessment

What the researcher cares about:

- Yield – did you get the number of reads you expected?
- Error - are the bases of reliable quality?
- Representative of sample – do the reads accurately represent the sample?

How do we get this information?

1. Think about sample quality and library quality

- Input sample quality is crucial for good sequencing
- High quality
- If quantity is low, use a kit designed for low input
- Check
 - spectrophotometric (Nanodrop)
 - fluorimetric (Pico- and Ribo-Green)
 - gel electrophoretic methods (Bioanalyser)
 - RNA Integrity Number (RIN)
- If you have someone else prepare your libraries, they should be able to tell you minimum necessary quantity, concentration, and QC values.



1. Think about sample quality and library quality

- Library quality

- Bioanalyzer
- Fragment Analyzer
- Agarose Gel

<https://support.illumina.com/bulletins/2016/05/library-quantification-and-quality-control-quick-reference-guide.html>

- Library quantity


- Fluorometric = a dsDNA-specific fluorescent dye method, such as QuBit, PicoGreen, and AccuClear
- qPCR

- Why?

- confirm insert size
- no primer dimers


2. Instrument Metrics

- Cluster density
 - You need the right amount of DNA per run, which varies by instrument
 - Under clustering – quality is generally fine but you lost yield
 - Over clustering – quality problems!
 - Avoid by properly quantifying your library




2. Instrument Metrics

- PF% - percent passing filtering
 - single molecule = single cluster = clear signal
 - Anything else doesn't pass the filter
- Phasing/Prephasing
 - Sometimes individual molecules in a cluster become out of sync
- Q30 – bases over quality value 30



PhiX Spike In

- Libraries must be diverse for proper sequencing
- Fix with a Phi-X spike in
 - A known quantity of DNA from a bacteriophage
 - 10-20% of sequences (or more)



Krueger et al 2011

<https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0016607>

3. Data Assessment

Software options

Illumina SAV – Sequence analysis viewer

- Part of BaseSpace
- Your sequencing facility probably uses this

FastQC

- Free
- Runs on linux



FastQC – Basic Statistics




Basic Statistics

Measure	Value
Filename	good_sequence_short.txt
File type	Conventional base calls
Encoding	Illumina 1.5
Total Sequences	250000
Sequences flagged as poor quality	0
Sequence length	40
%GC	45

FastQC – Per Base Sequence Quality

✓ Per base sequence quality


GOOD



FastQC – Per Base Sequence Quality

❌ **Per base sequence quality**

BAD



FastQC – Example Reports

- Good report

https://www.bioinformatics.babraham.ac.uk/projects/fastqc/good_sequence_short_fastqc.html

- Bad report


https://www.bioinformatics.babraham.ac.uk/projects/fastqc/bad_sequence_fastqc.html

- Adapter Dimer report (ligated adapters w/ no insert sequence)

https://www.bioinformatics.babraham.ac.uk/projects/fastqc/RNA-Seq_fastqc.html

- Small RNA with read through report:

https://www.bioinformatics.babraham.ac.uk/projects/fastqc/small_rna_fastqc.html




Trimming

Trimming


- From the quality control step, we know where the problems are
- All Illumina reads tend to have degrading quality at the end of the read
- Get rid of the bad data, keep the good data
 - Cut adapter sequences from the read.
 - Trim off low quality bases
 - Drop a read entirely if it is too low quality or too short

Adapter Trimming


Library Fragment:



Read returned from sequencing facility:




After Trimming:



Quality Trimming


cDNA
??

Quality Value




Quality Trimming

cDNA
??



What is your cut off?

Quality Trimming



Generally its ok to keep a single base or a few low quality bases to preserve downstream quality. But once the quality has degraded to the point of the bases being largely useless, make the trim.

Is trimming necessary?



- Depends on what you are doing with the data
- Balance data loss with downstream accuracy

If you meet all these criteria, maybe not:

- The reads are of high quality and have minimal adapter contamination
- You are mapping the reads to a well annotated reference genome
- You are doing gene quantification (no assembly, no variant calling)

Trimmomatic

- Optimized for Illumina NGS
- Very flexible
- Handles paired end data well
- Threaded
- Detects adapter read through
- No read through:
- Read through:



Bolger, A. M., Lohse, M., & Usadel, B. (2014). Trimmomatic: A flexible trimmer for Illumina Sequence Data. *Bioinformatics*, btu170.

Trimmomatic

Defaults are very stringent, but you can adjust.

The current trimming steps are:

- **ILLUMINACLIP**: Cut adapter and other illumina-specific sequences from the read. Must specify adapter sequence. **Comes with basic Illumina adapter files, make sure yours are in there or add yours!**
- **SLIDINGWINDOW**: Perform a sliding window trimming, cutting once the average quality within the window falls below a threshold.
- **LEADING**: Cut bases off the start of a read, if below a threshold quality
- **TRAILING**: Cut bases off the end of a read, if below a threshold quality
- **CROP**: Cut the read to a specified length
- **HEADCROP**: Cut the specified number of bases from the start of the read
- **MINLEN**: Drop the read if it is below a specified length after trimming

Skewer

- Faster than trimmomatic
- “Gentle” quality trimming by default
- Utilizes quality scores in adapter identification and allows insertions/deletions

<https://github.com/relipmoc/skewer>

Jiang, H., Lei, R., Ding, S.W. and Zhu, S. (2014) Skewer: a fast and accurate adapter trimmer for next-generation sequencing paired-end reads. BMC Bioinformatics, 15, 182.

Trimming

- Current community wisdom:
 - Quality trimming reduces error
 - But also reduces content and contiguity
- Gentle trimming is preferred – many times the defaults are too stringent, you will lose lots of data!
- Application matters
 - For expression, gentle to no trimming (phred 3 to 5)
 - For assembly and variant calling, trimming is good (phred 10 to 15)
 - Also read correction can make a difference!

More reading on trimming

- Williams et al. 2016 **Trimming of sequence reads alters RNA-Seq gene expression estimates**
- MacManes 2014 **On the optimal trimming of high-throughput mRNA sequence data**

More on read correction for transcriptome assemblies:

- Song and Florea 2015 Rcorrector: efficient and accurate error **correction** for Illumina **RNA**-seq reads
- MacManes and Eisen 2013 Improving transcriptome assembly through error correction of high-throughput sequence reads
- Heydari et al 2017 Evaluation of the impact of Illumina error correction tools on de novo genome assembly